

# COOL SOAK

## Stain & Rust Remover

Cool Soak™ Stain & Rust Remover is specially formulated to remove rust, water spots and detergent residue from stainless steel and tungsten carbide surgical instruments, automatic washing equipment and autoclaves.

### FEATURES & BENEFITS

- **Removes Stains and Rust.** Restores instruments' luster and shine.
- **Dissolves Residues.** Keeps box locks moving freely.
- **Instrument Safe.** If used properly, does not harm stainless steel or tungsten carbide instruments.
- **Bacteriostatic.** If stored properly, not favorable to bacterial growth.



**ORDER INFORMATION** ▶ Item Code: T5 • Size: 1 Gallon (3.8 Liters) • Case: 4 Units/Case

# COOL SOAK

Stain & Rust Remover



**USES** Manual • Ultrasonic • Automatic Washer • Cart Washer

**PHYSICAL PROPERTIES** Form: Liquid • Color: Violet • Odor: Mild • pH Concentration: <1

**TEMPERATURE** Usage: Room Temperature • Storage: Cool, dry area

## DIRECTIONS

Wear appropriate personal protective equipment based on safe work practices.

### MANUAL DESCALING

Dilute 1 part Cool Soak with 7 parts warm water. Soak thoroughly cleaned and opened instruments in prepared solution for 5-15 minutes. Scrub instruments under water with the Cool Soak solution and rinse thoroughly with water. Thoroughly clean instruments again once Cool Soak has been used and then lubricate and sterilize as normal. Change solution daily or more often as needed.

### ULTRASONIC DESCALING

Dilute 1 part Cool Soak with 7 parts water. Descaling

### AUTOMATIC WASHING EQUIPMENT

Dilute 1 fl. oz. (30 mL) to 3 fl. oz. (88 mL) Cool Soak per gallon (3.8 L) of water and follow manufacturer's guidelines for descaling.

### DESCALING AUTOCLAVES

Dilute 1 part Cool Soak with 1 part water. Apply solution to cool chamber walls and soak for 10-15 minutes. Wipe off (scrubbing may be necessary for heavier deposits) and rinse thoroughly with warm water.