

# OPTI-CIDE<sup>3</sup>

## Instrument Soak with Rust Inhibitor

Use this product as a disinfectant instrument soak. Opti-Cide<sup>3</sup>® effectively and quickly loosens and suspends blood and bodily fluids, non-organic and organic soil, and other cellular debris for easy removal. Disinfectant and cleaner for surgical instruments and other similar surfaces.

Opti-Cide<sup>3</sup> is on the [EPA's List N](#) indicating it can be used against SARS-CoV-2, the virus that causes COVID-19

### FEATURES & BENEFITS

- 2-Minute hard, non-porous surface disinfectant
- Bactericidal, Fungicidal, Virucidal, Tuberculocidal
- Ready-to-Use



**ORDER INFORMATION** ▶ Item Code: OCP04-128RI • Size: 1 Gallon (3.8 L) • Case: 4 Units/Case

## MEDICAL / SURGICAL / HOSPITAL USE

**2**  
MINUTE  
KILL TIME

Acinetobacter baumannii	Listeria monocytogenes
Avian Influenza A (H3N2) Reassortant Virus	Mycobacterium bovis, BCG (tuberculosis) (TB)
Escherichia coli - with ESBL (E. coli)	Pseudomonas aeruginosa
Hepatitis B Virus (HBV)*	Rotavirus Strain WA
Hepatitis C Virus (HCV)*	Salmonella enterica
Herpes Simplex Virus (HSV 2)	Serratia marcescens
Human Immunodeficiency Virus (HIV-1)*	Staphylococcus aureus
Swine Influenza A Virus (H1N1)	Methicillin-Resistant Staphylococcus aureus (MRSA)
Influenza A Virus Strain A2/Hong Kong	Trichophyton mentagrophytes (ringworm)
Carbapenem-resistant Klebsiella pneumoniae (CRKP)	Vancomycin-Resistant Enterococcus faecalis (VRE)

## VETERINARY / PET CARE FACILITY USE

**3**  
MINUTE  
KILL TIME

Bordetella bronchiseptica (kennel cough)	Chlamydia psittaci (parrot fever)
Bovine Viral Diarrhea Hepatitis C Virus	Feline Leukemia Virus
Duck Hepatitis B Virus	Infectious Bovine Rhinotracheitis Virus
Canine Distemper Virus	Murine Norovirus
Canine Para-Influenza Virus	Streptococcus equi

### HOSPITAL AND GENERAL DISINFECTION DIRECTIONS FOR USE

#### DIRECTIONS

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Apply this undiluted product by sponge, rag, cloth, mop, spray or immersion to wet all surfaces thoroughly.

#### FOR USE AS A PRE-CLEANING INSTRUMENT SOAK

Fill appropriately-sized container with a sufficient amount of this product to allow for complete submersion of items. Place items into solution and allow to soak for 2 minutes. Be certain to flush lumens of hollow devices with this product. Remove items and rinse. For critical and semi-critical items, follow with the appropriate cleaning and disinfection or sterilization process. Change solution as needed, i.e., when the solution becomes dirty or visibly soiled.

#### FOR USE AS A NON-CRITICAL MEDICAL ITEM DISINFECTION IMMERSION SOLUTION

Thoroughly pre-clean items being sure to remove all gross filth and heavy soil, then rinse and dry. Allowing for complete submersion of item(s), place pre-cleaned item(s) into an appropriate container that has been filled with this product. Flush lumens of hollow devices with this product. Allow item(s) to remain submerged for 2 minutes. Remove, rinse with clean, potable water and dry. Change solution daily or when it becomes visibly soiled or discolored.

### VETERINARY/PET CARE DIRECTIONS FOR USE

#### DIRECTIONS

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Apply this undiluted product by sponge, rag, cloth, mop, spray or immersion to wet all surfaces thoroughly.

#### FOR USE AS A PRE-CLEANING INSTRUMENT SOAK

Fill appropriately-sized container with a sufficient amount of this product to allow for complete submersion of items. Place items into solution and allow to soak for 3 minutes. Be certain to flush lumens of hollow devices with this product. Remove items and rinse. For critical and semi-critical items, follow with the appropriate cleaning and disinfection or sterilization process. Change solution as needed, i.e., when the solution becomes dirty or visibly soiled.

#### FOR USE AS A NON-CRITICAL MEDICAL ITEM DISINFECTION IMMERSION SOLUTION

Thoroughly pre-clean items being sure to remove all gross filth and heavy soil, then rinse and dry. Allowing for complete submersion of item(s), place pre-cleaned item(s) into an appropriate container that has been filled with this product. Flush lumens of hollow devices with this product. Allow item(s) to remain submerged for 3 minutes. Remove, rinse with clean, potable water and dry. Change solution daily or when it becomes visibly soiled or discolored.

This product cannot be used to disinfect critical or semi-critical medical devices.