

# OPTI-CIDE<sup>3</sup><sup>®</sup>

## Instrument Soak with Rust Inhibitor



Use this product as a disinfectant instrument soak. Opti-Cide<sup>3</sup><sup>®</sup> effectively and quickly loosens and suspends blood and bodily fluids, non-organic and organic soil, and other cellular debris for easy removal. Disinfectant and cleaner for surgical instruments and other similar surfaces.

- **2-Minute hard, non-porous surface disinfectant**
- **Tuberculocidal, Virucidal, Fungicidal, Bactericidal**
- **Ready-to-Use**

Order No.	Size
OCP04-128RI	1 Gallon (3.8 L)
OCPRI144-002	2 fl. oz. (59 mL)

## Instrument Soak with Rust Inhibitor

### Medical/Surgical/Hospital Use: 2-minute contact time

*Acinetobacter baumannii*  
 Avian Influenza A Reassortant Virus (H3N2)  
*Escherichia coli* – with ESBL (E. coli)  
 Hepatitis B Virus (HBV)  
 Hepatitis C Virus (HCV)  
 Herpes Simplex Virus (HSV-2)  
 Human Immunodeficiency Virus (HIV-1)  
 Swine Influenza A Virus (H1N1)  
 Influenza A Virus Strain A2/Hong Kong  
 Carbapenem-resistant *Klebsiella pneumoniae* (CRKP)  
*Listeria monocytogenes*  
*Mycobacterium bovis* BCG (tuberculosis) (TB)  
*Pseudomonas aeruginosa*  
 Rotavirus Strain WA  
*Salmonella enterica*  
*Serratia marcescens*  
*Staphylococcus aureus*  
 Methicillin-Resistant *Staphylococcus aureus* (MRSA)  
*Trichophyton interdigitale* (ringworm)  
 Vancomycin-Resistant *Enterococcus faecalis* (VRE)

### Veterinary/Pet Care Facility Use: 3-minute contact time

*Bordetella bronchiseptica* (kennel cough)  
 Bovine Viral Diarrhea Hepatitis C Virus  
 Duck Hepatitis B Virus  
 Canine Distemper Virus  
 Canine Para-Influenza Virus  
*Chlamydia psittaci* (parrot fever)  
 Feline Leukemia Virus  
 Infectious Bovine Rhinotracheitis Virus  
 Murine Norovirus  
*Streptococcus equi*

### Procedures Hospital and General Disinfection Directions for Use

<b>Directions:</b>	It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Apply this undiluted product by sponge, rag, cloth, mop, spray or immersion to wet all surfaces thoroughly.
<b>For use as a pre-cleaning instrument soak:</b>	Fill appropriately-sized container with a sufficient amount of this product to allow for complete submersion of items. Place items into solution and allow to soak for 2 minutes. Be certain to flush lumens of hollow devices with this product. Remove items and rinse. For critical and semi-critical items, follow with the appropriate cleaning and disinfection or sterilization process. Change solution as needed, i.e., when the solution becomes dirty or visibly soiled.
<b>For use as a non-critical medical item disinfection immersion solution:</b>	Thoroughly pre-clean items being sure to remove all gross filth and heavy soil, then rinse and dry. Allowing for complete submersion of item(s), place pre-cleaned item(s) into an appropriate container that has been filled with this product. Flush lumens of hollow devices with this product. Allow item(s) to remain submerged for 2 minutes. Remove, rinse with clean, potable water and dry. Change solution daily or when it becomes visibly soiled or discolored.

### Procedures Veterinary/Pet Care Directions for Use

<b>Directions:</b>	It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Apply this undiluted product by sponge, rag, cloth, mop, spray or immersion to wet all surfaces thoroughly.
<b>For use as a pre-cleaning instrument soak:</b>	Fill appropriately-sized container with a sufficient amount of this product to allow for complete submersion of items. Place items into solution and allow to soak for 3 minutes. Be certain to flush lumens of hollow devices with this product. Remove items and rinse. For critical and semi-critical items, follow with the appropriate cleaning and disinfection or sterilization process. Change solution as needed, i.e., when the solution becomes dirty or visibly soiled.
<b>For use as a non-critical medical item disinfection immersion solution:</b>	Thoroughly pre-clean items being sure to remove all gross filth and heavy soil, then rinse and dry. Allowing for complete submersion of item(s), place pre-cleaned item(s) into an appropriate container that has been filled with this product. Flush lumens of hollow devices with this product. Allow item(s) to remain submerged for 3 minutes. Remove, rinse with clean, potable water and dry. Change solution daily or when it becomes visibly soiled or discolored.

This product cannot be used to disinfect critical or semi-critical medical devices.