

Disinfectant Cleaner – Gallon Liquid

Opti-Cide[®] MAX cleans and disinfects hard, non-porous surfaces such as acrylic, metals (aluminum, brass, copper, stainless steel, chrome), glass, plastics (polycarbonate, polypropylene, polystyrene and polyvinylchloride (PVC)), Plexiglas[®], and vinyl. Opti-Cide[®] MAX kills hard-to-kill ESKAPE and ESKAPE MDRO organisms with a 1 minute ♦ kill time.

Opti-Cide MAX is on the [EPA's List N](#) indicating it can be used against SARS-CoV-2, the virus that causes COVID-19

FEATURES & BENEFITS

- **1 Minute ♦ Kill Time.** Kills tested organisms in 30 seconds to 1 minute. Faster kill times mean faster turnovers.
- **Kills ESKAPE and ESKAPE MDRO.** Kills dangerous ESKAPE and ESKAPE MDRO organisms in 1 minute. ESKAPE bacteria are responsible for approximately 2/3 of all HAIs.
- **Broad Spectrum.** A true broad spectrum disinfectant that is Bactericidal, Fungicidal, Virucidal, and Tuberculocidal.
- **Surface Safe.** Safe on clinical surfaces such as acrylic, aluminum, brass, copper, stainless steel, chrome, plastics and more. Maintain clear, bright plastics that are typically clouded by high levels of alcohol.
- **Low Alcohol Formula.** Cleans better than high alcohol formulas (55%). Low alcohol formula (20%) won't affix blood and soil to surfaces.



ORDER INFORMATION ▶ Item Code: M60035 • Size: 1 Gallon (3.8 L) • Case: 4 Units/Case

Disinfectant Cleaner – Gallon Liquid

DIRECTIONS

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

TO USE AS A CLEANER AND DEODORIZER

Apply solution to surface and allow surface to remain wet for 30 seconds prior to wiping with a paper towel or clean towel/sponge. This product can also be used to pre-clean all surfaces to be disinfected.

TO DISINFECT HARD, NON-POROUS SURFACES/ARTICLES

Spray solution to thoroughly wet pre-cleaned hard, non-porous surfaces. Allow surface to remain wet for 1 minute ♦ prior to wiping or allow surfaces to air dry.

Note: Disinfected surfaces that may come in contact with food must be rinsed with potable water after cleaning and/or disinfection and before reuse.

♦ This symbol indicates a 2 minute contact time is required to kill TB.

Kills ESKAPE Organisms

1 minute contact time

- [E] Enterococcus faecium
- [S] Staphylococcus aureus
- [K] Klebsiella pneumoniae
- [A] Acinetobacter baumannii
- [P] Pseudomonas aeruginosa
- [E] Enterobacter aerogenes

Kills ESKAPE MDRO Organisms

1 minute contact time

- [E] Enterococcus faecium
- [S] Methicillin-Resistant Staphylococcus aureus (MRSA)
- [K] Carbapenem Resistant Klebsiella pneumoniae NDM
- [A] Multi-Drug Resistant Acinetobacter baumannii
- [P] Pseudomonas aeruginosa
- [E] Enterobacter cloacae