

# OPTI-CIDE<sup>3</sup><sup>®</sup>

## Hard, Non-Porous Surface Sanitizer\*\* Disinfectant Cleaner



Use this product to clean, disinfect and sanitize\*\* hard, non-porous non-food contact, inanimate surfaces such as those made of acrylic, metals (aluminum, brass, stainless steel, chrome), glass, plastics (polycarbonate, polypropylene, polystyrene and polyvinylchloride), clear plastics, Plexiglas<sup>®</sup>, vinyl as well as laminated and painted surfaces. Opti-Cide<sup>3</sup><sup>®</sup> effectively and quickly loosens and suspends blood and bodily fluids, non-organic and organic soil, and other cellular debris for easy removal.

- **2-Minute hard, non-porous surface disinfectant**
- **10-Second soft surface sanitizer\*\* (Solution only)**
- **Bactericidal, Fungicidal, Virucidal\*, Tuberculocidal**
- **Ready-to-Use**

Order No.	Size	Order No.	Size
<b>OCS12-024</b>	24 fl. oz. (710 mL) Spray	<b>OCP01-1920</b>	15 Gallon (57 L)
<b>OCP12-024</b>	24 fl. oz. (710 mL) Pour	<b>OCP01-3840</b>	30 Gallon (114 L)
<b>OCP04-128</b>	1 Gallon (3.79 L)	<b>OCP01-7040</b>	55 Gallon (208 L)
<b>OCP02-320</b>	2.5 Gallon (9.5 L)	<b>OCW06-100</b>	100 Wipes

## Hard, Non-Porous Surface Sanitizer\*\* Disinfectant Cleaner

### Medical/Surgical/Hospital Use: 2-minute contact time

*Acinetobacter baumannii*  
Avian Influenza A (H3N2) Reassortant Virus  
*Escherichia coli* – with ESBL (E. coli)  
Hepatitis B Virus (HBV)\*  
Hepatitis C Virus (HCV)\*  
Herpes Simplex II Virus (HSV 2) G Strain  
Human Immunodeficiency Virus (HIV-1)\*  
Swine Influenza A Virus (H1N1)  
Influenza A Virus Strain A2/Hong Kong  
*Klebsiella pneumoniae*- Carbapenem-resistant (KPC)  
*Listeria monocytogenes*  
*Mycobacterium bovis*, BCG (tuberculosis) (TB)  
*Pseudomonas aeruginosa*  
Rotavirus Strain WA  
*Salmonella enterica*  
*Serratia marcescens*  
*Staphylococcus aureus*  
Methicillin-Resistant *Staphylococcus aureus* (MRSA)  
*Trichophyton mentagrophytes* (ringworm)  
Vancomycin-Resistant *Enterococcus faecium* (VRE)

### Veterinary/Pet Care Facility Use: 3-minute contact time

*Bordetella bronchiseptica* (kennel cough)  
Bovine Viral Diarrhea Virus  
Duck Hepatitis B Virus  
Canine Distemper Virus  
Canine Para-Influenza Virus  
*Chlamydia psittaci* (parrot fever)  
Feline Leukemia Virus  
Infectious Bovine Rhinotracheitis Virus  
Murine Norovirus  
*Streptococcus equi*  
Canine Parvovirus (5 minute contact time – wipes only)

Procedures	Solution Directions for Use
<p><b>Directions:</b></p> <p><b>NON-FOOD CONTACT SURFACES:</b></p> <p><b>For surface cleaning:</b></p> <p><b>For surface sanitization:</b></p> <p><b>For surface disinfection:</b></p>	<p>It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Apply this undiluted product by sponge, rag, cloth, mop, spray or immersion to wet all surfaces thoroughly.</p> <p>Prior to surface disinfection or sanitization**, pre-clean surfaces, making sure to remove all gross filth and heavy soil, by saturating surfaces with this product by spray, mop, rag, towel, paper towel, pour or immersion then wipe surface clean. Repeat, if necessary, until surfaces are visibly clean. Discard paper towels or launder cloth towels before reusing.</p> <p>Apply this product to clean surfaces by sponge, cloth, mop, spray or immersion to wet all surfaces thoroughly. Allow to remain visibly wet for 10 seconds for non-food contact surfaces at room temperature, then wipe surfaces using a clean mop, paper or cloth towel, or rinse surface using potable water or wipe surface dry or allow to air dry. Discard or launder spent materials before reuse.</p> <p>Apply this product to pre-cleaned surfaces by sponge, cloth, mop, spray or immersion to wet all surfaces thoroughly. Allow to remain wet for the appropriate time indicated for the purpose intended, then wipe dry using a mop, clean paper or cloth towel, or rinse surface using potable water and wipe surface dry or allow to air dry. Discard or launder spent materials before reuse.</p> <p>** When indicated for Sanitization use (for non-food contact surfaces)</p>

Procedures	Wipes Directions for Use
<p><b>Dispensing Directions</b></p> <p><b>Cleaning and Sanitizing/Disinfection of Surfaces:</b></p>	<p>To start feed, remove cover and discard seal. From the center of the wipe roll, pull up a wipe corner, twist it into a point and thread it through the hole located on the container cover. Pull through about one inch. Replace cover. Pull out first wipe and tear off at an angle. Remaining wipes feed automatically, ready for the next use. When not in use, keep center cap closed to prevent evaporation.</p> <p>Use a fresh wipe to pre-clean surfaces of all gross filth and heavy soil. Repeat as necessary until surfaces are visibly clean. To effectively sanitize or disinfect the pre-cleaned surfaces, use a fresh wipe or turn the wipe over to its clean side to thoroughly wet the surfaces and allow surfaces to remain wet for the appropriate time indicated for the purpose intended.</p>