

Product Name: Micro-Scientific Enzyclean II LS

Product Code: Z6

Ingredients	CAS No.	Function
Water	7732-18-5	Diluent
Glycerin	56-81-5	Stabilizer
Sodium hydroxide	1310-73-2	Buffer
Protease enzyme	9014-01-1	Cleaning agent
Amylase enzyme	9000-90-2	Cleaning agent
Methylisothiazolinone	2682-20-4	Preservative
Benzisothiazolinone	2634-33-5	Preservative
Acid blue 9	3844-45-9	Colorant
Yellow 5	1934-21-0	Colorant
Octoxynol	9036-19-5	Fragrance
d-Limonene	5989-27-5	Fragrance