

**Product Name:** Opti-Cide Max Wipes Disinfectant Cleaner

**Product Code:** M60034, M60038

<b>Ingredients</b>	<b>CAS No.</b>	<b>Function</b>
Water	7732-18-5	Diluent
Isopropanol	67-63-0	Antimicrobial
Ethanol	64-17-5	Antimicrobial
Butoxyethanol	111-76-2	Solvent
Alkyl C12-18 dimethylbenzyl ammonium chloride	68391-01-5	Antimicrobial
Alkyl C12-14 dimethyl ethylbenzyl ammonium chloride	85409-23-0	Antimicrobial
Propylene glycol	57-55-6	Solvent
Tetrasodium EDTA	64-02-8	Chelating agent
Amine oxide	Withheld	Cleaning agent
Benzotriazole	95-14-7	Corrosion inhibitor
Inorganic base	Withheld	pH adjuster
Amine oxide	Withheld	Cleaning agent