

Product Name: Micro-Scientific Enzyclean II NS

Product Code: ENF1, ENF5, ENF32

Ingredients	CAS No.	Function
Water	7732-18-5	Diluent
Glycerin	56-81-5	Stabilizer
Alkanoate salt	764-71-6	Cleaning agent
Sodium hydroxide	1310-73-2	Buffer
Amylase enzyme	9000-90-2	Cleaning agent
Protease enzyme	9014-01-1	Cleaning agent
Methylisothiazolinone	2682-20-4	Preservative
Benzisothiazolinone	2634-33-5	Preservative
Acid blue 9	3844-45-9	Colorant
Yellow 5	1932-21-0	Colorant
Octoxynol	9036-19-5	Fragrance
d-Limonene	5989-27-5	Fragrance