

**Product Name:** Opti-Cide<sup>3</sup>

**Product Code:** OCP01-1920, OCP01-3840, OCP01-7040, OCP02-320, OCP04-128, OCS12-008, OCS12-024, S5-1, S5Q-1

<b>Ingredients</b>	<b>CAS No.</b>	<b>Function</b>
Water	7732-18-5	Diluent
Isopropanol	67-63-0	Antimicrobial
Butoxyethanol	111-76-2	Solvent
Alkyl C12-18 dimethylbenzyl ammonium chloride	68391-01-5	Antimicrobial
Alkyl C12-14 dimethyl ethylbenzyl ammonium chloride	85409-23-0	Antimicrobial
Ethanol	64-17-5	Solvent
Sulfate of ethoxylated alcohol	Withheld	Cleaning agent
Nonionic surfactant	Withheld	Cleaning agent